

Water Quality Test



BONDURANT
TECHNOLOGIES INC.

Date : 09/09/2020

Enzyme Substrate Coliform Test

Method:	Standard Method – SM9223-B (colilert-18)
Analysis Start:	15:08 09/09/2020
Analysis End:	09:08 09/10/2020
Incubator Start Temp.	34.8°F
Incubator End Temp.	35.1°F

Positive Wells Total Coliform

Pre-Filter:	Large - 49	Small - 48	>2,419.6 MPN/100mL
Post-Filter:	Large – 0 UNDETECTED	Small – 0 UNDETECTED	<1 MPN/100mL

Positive Wells E.coli

Pre-Filter:	Large - 49	Small - 48	>2,419.6 MPN/100mL
Post-Filter:	Large – 0 UNDETECTED	Small – 0 UNDETECTED	<1 MPN/100mL